


# MODAS Instrumental Control and Data Analysis Software


## for MT Explorer 100 Ion Trap Mass Spectrometer

MODAS provides powerful and flexible instrumental control and data analysis for MassTech's MT Explorer 100 ion trap mass spectrometer. Supported operating systems are Windows Server 2003; Windows Server 2008; Windows 7, Windows Vista; and Windows XP. The software has 32-bit and 64-bit versions. MODAS has flexible modular architecture designed for easy maintainability and extensibility.




### Main features

- MS and MS/MS modes
- Normal and Expert mode of operations
- Multiple parent ions
- Automatic gain control
- Built-in ion source control modules
- GC and LC ready
- Real-time spectrum analysis
- Built-in chromatogram viewer



Normal and expert mode of operations



Built-in chromatogram viewer



6992 Columbia Gateway Drive  
Columbia, MD 21046, USA  
Telephone: 443.539.1758  
Fax: 443.539.1759  
Web: [www.apmaldi.com](http://www.apmaldi.com)  
Email: [msms@apmaldi.com](mailto:msms@apmaldi.com)