



MALDI + Ion Funnel = Sub-Atmospheric Pressure MALDI

SubAP/MALDI(ng)


A powerful tool for high-resolution MS imaging



- Operates over 1-10 Torr pressure range
- Variable laser spot size (10 to 50 micron)
- Easy exchangeable with the ESI source
- Available for ThermoFisher Q Exactive, LTQ Orbitrap



Other Ion Sources and Products



Ion Funnel Interface with exchangeable ESI/AP-MALDI sources



AP-MALDI (ng) UHR



MT Explorer 50
Field-Deployable Tandem MS with
DS-sESI and DS-APCI Ion Sources



MT Explorer 50
Field-Deployable Tandem MS
with DART Ion Source